MEBER<sup>+</sup>

# ERGON

Patent  
pending

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**en** Technology and advanced rescue techniques

[www.meber.it](http://www.meber.it)  
[www.meber.it/ergon](http://www.meber.it/ergon)

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Correct treatment of spinal injuries starts on the accident  
scene

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
## **Variable geometry comfort scoop stretcher**

Designed by rescue professionals and technicians with expertise in medical emergency devices for recovery in critical situations.

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New operative scenarios in emergency rescue

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### en Introduction

Spinal injuries are traumatic events that cause damage that often results in permanent invalidity. This kind of injury weighs heavily on public health with a high social and economic cost for modern society, as youngsters are the most frequently affected age group. The aim of in situ patient immobilization and the subsequent need for mobility and accessibility to the patient during pre-hospital rescue phases is to protect the spine from passive movements that could worsen the injury caused by the accident. In all seriously injured patients, spinal column involvement must be suspected until proven to the contrary.



Fig. 1

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

#### **en New perimetric protections**

Hard plastic protections are now present on ergon corners to avoid possible soft high comfort material damages during rescue operations and stowage. Moreover these protection grant a better sliding also inside rescue vehicles stowage compartments.

#### **New high visibility indications**

All indications how to operate the flexilock locking device are now available with colour contrast so to make them clear even in case of low visibility and lightning.

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**en continued from page 6**

One of the problems of emergency rescue is the increase in obese and grossly overweight patients. An ever-greater number of rescue patients are above standard weight or "bariatric patients". This is a significant problem for both rescue professionals and the equipment used, which have to be able to support great stress for long periods of time.

Domestic accidents are another increasingly frequent occurrence. This type of accident has reached a truly alarming rate both in this country and on an international scale. Accidents in the home, which interest all age brackets, are the first cause of death amongst children in developed countries, although housewives represent the most affected category overall. The greatest difficulty in this type of context is speed: the accident victim is often alone at home and unable to communicate with the emergency rescue switchboard immediately. This calls for versatile, easy-to-handle equipment. Unfortunately, accidents in the workplace make the news on a daily basis, due to the increase in the number of companies at environmental risk or risk of significant accidents.

Road accidents are the ninth cause of death worldwide and also in Italy they represent an important public health problem, as the first cause of death amongst the under-40s. The impact in terms of admissions to Accident and Emergency Departments and invalidity is considerable and it is estimated that every year this type of event leads to about one and a half million A&E admissions and about 20 thousand cases of invalidity.

Disaster scenario organisation and planning is a well-known issue and with the publication by the **Prime Minister's Cabinet Civil Protection Department** (Ministerial Decree 13th February 2003 - Adoption of the "Guidelines for medical rescue organisation in disasters") it has become a matter of great current interest and a demanding topic from an organisational, logistic and training standpoint.

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
Stretchers through history:  
the evolution of the species

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
1900

The first wood and fabric patient immobilization device was created by Neil Robertson in the early 1900s, it was of Japanese inspiration, was made of bamboo and stitched by hand.

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
1930

Designed in the 1930s, the Thomas stretcher is considered a landmark in mountain rescue. During the same period, the Mariner stretcher was also used by the Navy.

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
1950

For the British Everest Expedition in the 1950s, Donald Duff designed a steel tube stretcher.

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With the War in Vietnam, US army rescue crews started using transport devices that can be considered the forefathers of modern-day field stretchers. In the following years, Californian volunteer Rick Kendrick designed the KED stretcher for the extrication of accident victims.

From the 1980s and especially after the Irpinia earthquake, the scoop stretcher, considered by many to be the first real trauma device, was introduced on ambulances in Tuscany, Liguria and northern Italy in general.

The spinal board did not appear until after 1990. More recently, spinal boards have been designed with different sizes and using different materials and the KED has been restyled many times, although the basic concept and system of use have remained the same.

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
A new product for a new way  
of doing rescue

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## **en** The ideal device for rescue professionals

Given the evolution in the type of trauma and the experience gained in operating contexts, the stretchers of the future must guarantee practicality, availability, efficiency and rapidity during use and transportation. Ergon, Me.Ber.'s latest generation stretcher, has been designed to overcome the technical limits of similar devices present on the market. Specifically, it has been designed to simplify the work of rescue professionals, guarantee patient safety and comfort, operate in difficult situations, facilitate operator training, comply with current regulations and anticipate those of the future.

BARELLA ATRAUMATICA A  
GEOMETRIA VARIABILE  
VARIABLE GEOMETRY  
ATRAUMATIC STRETCHER

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## **en** The creation of a new product

In recent years, there has been a significant reduction in the use of scoop stretchers in favour of spinal boards and they are consequently now primarily used for rescue and very short transfers to other devices such as the spinal board or vacuum mattress. Ergon has been designed as a device that combines the ductility and flexibility of use of the scoop stretcher with the other characteristics of the spinal board.

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
Variable geometry  
comfort scoop  
stretcher

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Widespread stiffness with stress on points of contact (occipital area, shoulders, pelvis, calves and heels)

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The concave shape may restrict chest expansion. Made of metal and may not prevent hypothermia.

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Good level of comfort although it does not provide areas of differentiated softness.

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
The shape facilitates chest expansion. Comfort areas at main decubitus points.

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
Best device for spinal support and immobilization.

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Spinal support is precarious and incomplete.

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

Good spinal support. Minimal distance between the two halves of the stretcher.

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Possible access difficulties  
due to minimal distance from  
the ground.

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Potential interference of  
locking device with patient's  
head.

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Excellent handle distribution,  
all of the easy-grip type  
Off-centre locking devices.  
Handle area for longitudinal  
axis control.

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
**en** The greatest product for  
the smallest patients

Ergon puts its advanced technology at the service of young patients. The Junior version shares all the great benefits of the adult stretcher. The Flexilock locking device is off-centre to allow operation without invading the space occupied by the patient's head and the wider than standard blade delivers absolute comfort even to the most demanding patients.

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# IN DETTAGLIO IN DEPTH

ADVANCED TECHNOLOGY  
INNOVATIVE MATERIALS

INNOVATIVE VARIABLE  
GEOMETRY LOCKING DEVICE

Patent pending

DESIGNED FOR COMFORT,  
SAFETY AND ERGONOMICS

## **en** Heading for “best practice”

Ergon has set itself the mission of becoming the **best practice**, the landmark device in the emergency sector. The development and construction of Ergon incorporates some of the most advanced technologies ever used on a stretcher. Revolutionary shape, latest generation co-moulded materials and the innovative Flexilock locking device are just a few of the features that

make this product a candidate for excellence. They make Ergon an extremely comfortable stretcher for patients and safe and ergonomic for rescue professionals. Ergon's design makes it possible to restrict or eliminate risks for the patient deriving from long periods on the stretcher or as a consequence of inadequate handling.

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
The materials of tomorrow make  
the difference in emergency rescue

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**HC<sup>3</sup>** High Comfort Composite Cell

**en** Ergon is the first stretcher on the market in which materials with different mechanical characteristics are combined, using innovative processes. The aim is to create a more structurally resistant device.

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

Section of the layered  
coupling between plastic  
materials and metals.

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Section of  
high-resistance  
telescopic tube

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Section of coupling  
between locking device  
materials

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HIGH INSULATING  
THERMAL BARRIER

**en** Ergon is made using the injection moulding technique, with compound materials and metal inserts. It is composed of a soft material co-moulded onto a high resistance plastic polymer. The new system composed

of three special materials makes the stretcher resistant to collisions and changes of temperature and flame retardant (EN 1021).

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

X-ray and scan picture  
interpretation in  
standard projections  
in emergency  
conditions.

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The shape of the frame  
is slight and regular and  
does not cause confusion  
when interpreting  
radiographic pictures.

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

Images obtained with RM systems at 1,5T

## **en Urgent procedures in absolute safety**

In first aid it is often important that patient is subjected to radiological examination without being moved from the stretcher on which he arrived to first aid station: Ergon has been designed to allow doctors and

radiologists to evaluate with precision results of such tests. Device suitable to support patient during execution of radiological examinations (RX) in first aid field.

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The first variable geometry  
locking device

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Detail of the standard  
locking device in its new  
position and released  
variable geometry  
stretcher.

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Detail of the release  
device on the  
operating arm

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### **en** Different situations, a single solution

Ergon is fitted with a new generation locking device positioned off-centre and in opposite pairs along the two sides of the stretcher, to allow a variable geometry. This innovation allows the rescue professional to operate in most rescue situations made difficult due to the characteristics of the ground or to perform rescues in confined or hostile environments.

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All the flexibility needed for the most demanding rescue work

**"...PER RECUPERI IN AMBIENTI OSTILI..."**

"...FOR RESCUES IN HOSTILE ENVIRONMENTS..."

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**en** In situations in which the characteristics of the ground are difficult or handling space is restricted, use of a standard scoop stretcher is unsuitable, unless the patient is moved, which increases the risks of the lesions caused. In these cases, it is essential to use a device suitable for critical conditions.

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Detail of standard  
locking device opening

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
**en** Ergon is the only stretcher able to perform variable geometry manoeuvres. This particular configuration can be obtained using the releasing jointed arm of the Flexilock locking device. Ergon is easy to close even in difficult situations where conventional scoop stretchers are

unsuitable. Once the two half-stretchers have been closed, even if in a log-roll position, the patient rests solidly on the full length of the stretcher. The patient can now be transferred in safety and the rescue professionals can complete the rotation manoeuvre (complete supination).

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Anyone who knows how to use a scoop stretcher can use Ergon

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## **en** The same, but different

In standard situations, Ergon is used in exactly the same way as any known type of stretcher used for mobilisation: this helps even less-expert rescue professionals because the opening locking device is easy to identify and the stretcher is used as normal, without needing any particular skills or specific training.

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VERY HIGH TECHNOLOGY  
INNOVATIVE MATERIALS

**GANCIO INNOVATIVO A  
GEOMETRIA VARIABILE**

INNOVATIVE VARIABLE GEOMETRY  
LOCKING DEVICE

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SHAPES DESIGNED FOR COMFORT,  
SAFETY AND ERGONOMICS

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**en** The Flexilock locking device is positioned off centre to allow rescue professionals to perform the closure manoeuvre away from the areas occupied by the patient's head and feet and also to close the device, thus feeling the patient's weight to a lesser degree as it is positioned centrally. The Flexilock locking device is counterlateral to avoid accidental opening of the device. The Flexilock locking device's central non-slip manoeuvre handle makes it possible to move the stretcher sideways or allows the loading of heavy patients by a team of rescue professionals.

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Separation manoeuvre identical to that  
using a scoop stretcher

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High heat retention and insulation

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Patient's shoulder line inside the stretcher handling area

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Thermal insulation for rescue professionals who, thanks to the HC<sup>3</sup> handles, are able to work wearing just latex gloves, even at low temperatures.

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
Thanks to the revolutionary  
material HC<sup>3</sup>, Ergon  
guarantees maximum  
comfort even for long  
periods

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SHAPE OF THE BLADE  
FOR MAXIMUM  
RESTRAINT

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The special shape of the blade  
facilitates insertion underneath  
the patient

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
Maximum thermal insulation for  
the patient, obtained using the  
innovative HC<sup>3</sup> material

**en** The shape of the stretcher favours patient accommodation, accompanies the natural curve of the ribcage and favours restraint and correct respiratory compliance. The use of HC<sup>3</sup>, in addition to guaranteeing absolute structural stability, also allows high thermal

insulation. This helps prevent hypothermia and, thanks to the stretcher's shape, makes it possible to prolong the time the patient spends on the device without further stress.

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
Shape of the blade optimised for maximum comfort

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**Sagoma barella**  
Shape of the stretcher

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► Pag. 56

OPERATING PROCEDURES  
RESCUE WITH ERGON

**en** New technologies have allowed us to drastically reduce the thickness of the blade, guaranteeing adequate support and better resistance over a vaster contact surface. The over-moulding between the rigid


plastic module and the aluminium tube with an elliptic section, creates a single body that guarantees resistance to the various types of dynamic and static stress.

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SHAPES DESIGNED  
TO EXCEED THE LIMITS

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
**en** The blade is just 5 mm thick for facilitated insertion between the patient and the ground. It has an elongated form to guarantee better support of the spinal column. This shape, which has been tested on the field, prevents unintentional gathering of materials (grass, clothes, sheets, blankets, etc.) underneath the patient during loading.

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
Maximum resistance for the toughest of tasks

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## en Regulation-compliant

The tests performed on Ergon satisfy EN 1865 standards as regards load, flexion and torsion stress.

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With a patient weighing  
100kg, the gap between  
the two half stretchers  
increases by just 10 mm

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### **en** Bariatric rescue


Ergon satisfies the increasingly frequent demands for devices suitable for rescuing patients who exceed standard weight.

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
**en** Ergon is designed to aid rescue professionals with unique features. Around the perimeter frame, which is composed of tubes with an elliptic section, the handle grips are coated in HC<sup>3</sup> and rotated 20° to favour a better grasp. Special supports and eyelets have been positioned for easy fastening of belts with spring catches. There are also a number of slots for patient immobilization using conventional devices.

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Easy grip handles for absolute safety

Preso sul gancio  
Flexilock in HC<sup>3</sup>

HC<sup>3</sup> handle on Flexilock  
locking device

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Easy grip side handles  
for a safe grasp

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Non-slip handle for lifting  
off the ground

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**en** Ergon guarantees a multitude of grip positions, to suit requirements, the space available, the patients' size and weight, rescue crew habits and varying number of rescue professionals.


# ESTENSIONE MODULO PIEDI A SGANCIO RAPIDO

Extension of the  
quick-release footplate

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**en** The high visibility release devices identify the foot side and allow elongation to suit the patient's height. The stretcher reaches a length of 201 cm when fully extended.

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
High-visibility  
release devices

# X-LOCK, FERMACAPO UNIVERSALE ULTRA-PORTATILE

X-lock, ultra-portable universal headlock

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**en** Device featuring an innovative design composed of a single element with an integrated fastening belt to allow rapid fastening of the head to the stretcher, without further manipulation of the cervical spine. Made of multiply polyester and varying density expanded polyurethane, the X-lock headlock guarantees rapidity and optimisation of equipment/crew resources, speed and efficacy in decontamination and limitation of the "restraint discomfort" effect in the immobilised patient.

# X-DOUBLE, FERMACAPO ADATTABILE

X-Double, Padded adjustable head immobilizer

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**en** X-DOUBLE is a divisible and adjustable head immobilizer designed in two pieces that can be fixed to the stretcher with special straps. This kind of system allows to fix the head immobilizer to the stretcher and to separate the two parts of the stretcher in order to load the patient. Made of expanded material with closed cells, an inside aluminium sheet and covered by a liquid-proof and bacterial -proof product. Designed to be used also in presence of a motorbike or automotive safety helmet.

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
X-fix, belt restraint system with 10 attachment points

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**en** Standard restraint belts fitted with spring catches: innovative system designed to speed up and facilitate the work of rescue professionals that maintains high comfort whilst guaranteeing spinal column stability.

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Reinforce along the perimeter made  
of coupled canvas

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EN566 ribbon ring for  
winch operations

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Single shoulder straps  
on each sides.

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**en** Designed to enable easy transport of the Ergon on the accident scene, this transport bag ensures a full safety usage even in the most critical rescue situation. It is equipped with 4 shoulder-straps for a safe grip on a single shoulder. The dismantled stretcher is stored inside along with all its accessories, everything fit in extremely reduced dimensions. Centrally positioned there is one ribbon ring certified EN565 which enables all type of winch operations.

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**en** In just 9 cm thickness, this bag can contain Ergon stretcher and all its accessories. It is designed to be used in warehouses, also hanged as a coat hanger, or to be used on board of rescue vehicles. Once the bag is open, it becomes a perfect ground-isolated (from water, mud and dust) and useful washable support for the preparation of medicines or to operate in safety with electro-medical devices.

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## **en** Ergon in action

The reconstructions shown here are intended to indicate an operating sequence in those situations in which the patient's critical position exalts the qualities of Ergon illustrated previously, but should not be considered as an alternative to standard mobilisation-immobilization procedures.

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
**INTERVENTO**  
INTERVENTION

### **en The right choice**

Accidents in the workplace with serious consequences are increasingly frequent. This calls for correct preparation and suitable operative instruments. If the accident victim is found lying with his/her back against a wall or any other structure preventing access to the spinal area,

it is important to perform rescue without moving the patient until he/she has been immobilised. This is why Ergon, with the revolutionary Flexilock locking device, is the right choice. No other device makes rescue in these conditions possible.

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**SVINCOLO**  
SEPARATION

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### **en Immobilization first**

It is not always easy to formulate a correct diagnosis of the injury, especially when the patient is unconscious or in a state of confusion. The rescue worker must know what device to use and the best way to use it. Trauma patient immobilization is a key phase of rescue and is an

operation that is often not straightforward or easy. The patient's height and weight and any breathing difficulties must be taken into consideration. This is why Ergon, with its special shape, offers maximum patient comfort even in the case of long periods on the device.

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Trauma patient in critical position

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
**INTERVENTO**  
INTERVENTION

### **en One step higher**


Another case of rescue in prohibitive conditions occurs when the trauma patient is lying with his/her spine astride the curb or on a step. In this case, the patient must not be moved until the cervical spine has been completely immobilised. Moving the patient before

the operation has been concluded could worsen spinal damage. In such cases, a conventional scoop stretcher is unsuitable because it is not possible to guarantee patient immobility during the half-stretchers coupling phase.

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**SVINCOLO**  
SEPARATION

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### **en** When time means everything

When rescuing patients with suspected cervical injuries, the time factor is paramount to correct diagnosis and relative therapy. The use of a conventional scoop stretcher and being unable to move the patient until complete immobilization, is likely to prolong the time

needed for rescue and transfer to hospital. This could hinder the efforts of rescue workers and doctors. With Ergon, this kind of rescue takes less than three minutes, during which time the patient does not need to be moved.

# APPENDICE

## Appendix

SOCCORSO E DIRITTO (REFERRED TO ITALIAN REGULATIONS)

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
## CUSTOMISATION

Ergon is ready to be customised with your logo


### **it** Me.Ber quality with your name on it

Ergon features a number of raised areas specifically intended for identifying the owner of the device (stretchers often get switched between the various hospitals or ambulance associations by mistake or due to carelessness). Subject to request, the name of sponsors, patrons or device donors can be included.

# IL SISTEMA ERGON. GAMMA PRODOTTI.


The Ergon system. Product range.

## Art.: 16150/V

Barella comfort a cucchiaino a geometria variabile verde/blu

Variable geometry comfort scoop stretcher green/blue

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
Variable geometry comfort scoop stretcher yellow/blue

## Art.: 16160/V

Ergon Junior, barella comfort a cucchiaino a geometria variabile verde/blu

Ergon Junior, variable geometry comfort scoop stretcher green/blue

## Art.: 16160/G

Ergon Junior, barella comfort a cucchiaino a geometria variabile giallo/blu

Ergon Junior, variable geometry comfort scoop stretcher yellow/blue

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

# Index


<b>1</b>	<b>Introduction</b>	<b>6</b>
	New operative scenarios in emergency rescue	6
	Stretchers through history: the evolution of the species	10
	A new product for a new way of doing rescue	12
<b>2</b>	<b>Comparing technologies</b>	<b>14</b>
	Patient comfort	16
	Spinal support	17
	Rescue crew ergonomics	18
	Ergon junior	19
<b>3</b>	<b>Ergon in depth</b>	<b>20</b>
	The materials of tomorrow make the difference in emergency rescue	22
	X-ray compatibility properties	26
	Flexilock	28
	All the flexibility needed for the most demanding rescue work	30
	Variable geometry	32
	Anyone who knows how to use a scoop stretcher can use Ergon	34
	High heat retention and insulation	38
	Shape of the blade optimised for maximum comfort	40
	Maximum resistance for the toughest of tasks	44
	Maximum rescue crew ergonomics	46
	Easy grip handles for absolute safety	48
	Extension of the quick-release footplate	50
	X-lock, ultra-portable universal headlock	52
	X-double, padded adjustable head immobilizer	53
	X-fix, belt restraint system with 10 attachment points	54
	Borsa zainabile per il trasporto	56
<b>4</b>	<b>Ergon on the field</b>	<b>58</b>
	Occupational accident	60
	Trauma patient in critical position	64
<b>5</b>	<b>Appendix</b>	<b>68</b>
	The Ergon system. Product range.	70

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All rights reserved. Variations can be done without notice. Me.Ber. S.r.l. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.