



TECHNICAL DATASHEET

Item code	679	
Name	"EXTRA LITE ZENA" Reduced width chair, 4 foldable wheels with footrest, armrest and aluminium platform	
Dimensions [mm]	Min length	650
	Max length	1220
	Width	470 (rear levers closed) 500 (rear levers open)
	Height	1024
	Folded thickness	290
	Seat height	560
Used materials	Al alloys	
Wheels	n°2 Ø200 mm (fixed) n°2 Ø125 mm (rotating)	
Weight [kg]	13	
Max load [kg]	230	
Declared use	Manual transportation of a seated patient with use-compatible pathologies who are temporarily or permanently disabled	
Conformity	MDD 93/42/CEE – Class I Medical Device	



2007

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it