



## SCHEDA TECNICA

<b>Item code</b>	<b>7097/4RG PROOF</b>	
<b>Name</b>	<b>"MERCURY LITE 5" Aluminium self loading stretcher, 5 positions with 4 rotating wheels EN 1865 certified – High version</b>	
<b>Dimensions [mm]</b>	Length	1970
	Length with folded trolley	1620
	Width	560
	Height	1000 – 900 – 810 – 685 – 390
	Folded height	310
	Loading height	680 ÷ 700
<b>Used materials</b>	Aluminium alloys	
<b>Wheels</b>	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm	
<b>Weight [kg]</b>	42	
<b>Max load [kg]</b>	220	
<b>Declared use</b>	Transporting, in the lying down position, of temporarily or permanently disabled patients with pathologies that allow this. The stretcher is suitable for transporting patients on vehicles appropriately prepared with fastening systems	
<b>Conformity</b>	MDR2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified	



2109

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: [info@meber.it](mailto:info@meber.it)