

TECHNICAL DATA SHEET



Item code	894/SNOW
Name	<i>"SNAKE SNOW" Independent chambers vacuum mattress with anti-slip system for snow</i>
Dimensions [mm]	Length: 2100 Width: <i>head zone</i> 800 <i>middle zone</i> 1000 <i>foot zone</i> 800
Used materials	Polyurethane
N° of handles	8
N° of patient restraint straps	4
Range temperature [°C]	-30 ÷ +70
Peso [kg]	7,5
Declared use	Temporary immobilisation of the body or for compatible pathologies and eventual transfer onto suitable devices for movement
Conformity	MDD 93/42/EEC Class I Medical Device