




## 9042 Galaxy-P

Orange spinal board with pins



### FEATURES

Same features as “Galaxy” spinal boards Art. 9010, 9012 and 9015 equipped with 10 plastic PINS useful for a fast fixing of restraint systems with karabiners. Supplied with 3 blets (605/MEB).

### MAIN FEATURES

Length (cm)	183
Width (cm)	41
Thickness (cm)	5
Loading capacity (kg)	150
Weight (kg)	5,8
Trauma	<input checked="" type="checkbox"/>
Usage	Immobilization
Material	ABS

### ACCESSORIES



**Jelly 629**

Orange standard head immobilizer  
for spinal board