




9050 Nilix


Yellow basic spinal board with pins



FEATURES

Basic spinal board made of high density plastic material, equipped with 8 PINS.
It's supplied with 3 belts.

MAIN FEATURES

Length (cm)	183
Width (cm)	40
Thickness (cm)	5
Loading capacity (kg)	150
Weight (kg)	6,5
Trauma	
Usage	Immobilization
Material	HDPE (High-density polyethylene)

ACCESSORIES



Jelly 629

Orange standard head immobilizer
for spinal board
