



# 1345

Empty emergency bag



## FEATURES

**An emergency bag made of Nylon 840.** It is compact but roomy and, thanks to its internal modularity, it allows personalized organization of the available space so that everything required for an emergency can be housed in the best way. It is made up of 4 identical pockets, 2 side ones, a front one and a back one, which are prepared and provided with supports and elastics for holding various medical material. It also has another small removable bag with a transparent front that allows the contents to be seen. Thanks to its construction, the linings can be completely removed to allow internal cleaning and disinfection. Can be carried by hand, over the shoulder or as a backpack. **Dimensions (LxHxD): 45 x 33 x 34 cm.** Colour: red.

## MAIN FEATURES

Width (cm)	45
Height (cm)	33
Depth (cm.)	34
Weight (kg)	1,8
Material	Nylon

---

IMAGES

---



---

SEE ALSO

---



**Rescuemed 1343**  
Rescue knapsack



**1365**  
Rescue bag Nevada



**1342**  
Emergency empty "Baby Bag"

---