



## 9236 Artis

Cervical collar set 6 pieces with carrying case



### FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.

### MAIN FEATURES

Trauma



### SEE ALSO



**Artis 9220**  
Rigid cervical collar in 1 piece size  
INFANT



**Artis 9222**  
Rigid cervical collar in 1 piece size  
PEDIATRIC



**Artis 9224**  
Rigid cervical collar in 1 piece size  
NECKLESS



**Artis 9226**  
Rigid cervical collar in 1 piece size  
SHORT



**Artis 9228**  
Rigid cervical collar in 1 piece size  
REGULAR



**Artis 9230**  
Rigid cervical collar in 1 piece size  
TALL