



## 1354

Red knapsack "115" with side pockets and compartments



### FEATURES

Knapsack made of red anti-tear nylon, double seams, reflective strips, reinforced shoulder straps and strong 10 mm zips. Supplied with front opening wide space inside of which are 4 rack pockets to rationally store various devices. Possibility to fix 2 litres oxygen bottle. Two wide lateral sections contain four bags in different colours allowing fast identification of various kits. PVC Water proof bottom.

**Dimensions (LxHxD): 44x56x26 cm**

### MAIN FEATURES

Width (cm)	44
Height (cm)	56
Depth (cm.)	26
Weight (kg)	4
Material	Nylon

### IMAGES



---

SEE ALSO

---



**1362**  
"Ambumed" rescue knapsack



**1366**  
Red knapsack "115" with side  
pockets and compartments  
COMPLETE



**Rescuemed 1343**  
Rescue knapsack

---

ACCESSORIES

---



**1354-001**  
Carrying bag for knapsack "115"

---