



## 11064

Manual suction unit



### FEATURES

#### HAND SUCTION UNIT

Designed to provide quick and effective pharyngeal and tracheal suction for emergency use.

The adjustable stroke knob provides different suction pressure.

Supplied with 250 ml bottle and 2 non sterile adult cannulas 17 cm (Ø 6 mm) and 19 cm (Ø 10 mm).

Max vacuum: -400 mmHg

Peak free airflow: >20 l/min

External diameter of catheter connection: 10 and 6 mm

Weight: 260 g

Size: 230x60x190 mm

### MAIN FEATURES

|             |      |
|-------------|------|
| Length (cm) | 23   |
| Height (cm) | 19   |
| Width (cm)  | 6    |
| Weight (kg) | 0,26 |

---

ACCESSORY OF

---



**1366**

Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---