



7018/A

HI-COM: Orange mattress made with expanded polyurethane at slow memory and fluid therapeutic gel



FEATURES

Fluid therapeutic gel mattress with 3 different sections of expanded polyurethane (one of which at slow memory). Outer raised sides to contain patient and two rigid zones head/ foot. Sacral/dorsal zone is made by two fluid therapeutic gel sheets (cm 40x40xH1,5) to create a balanced comfort for the patient. 3 expanded polyurethane zones absorb vibrations due to passive solicitation movement of vehicle, while two zones in fluid therapeutic gel absorb patient's weight granting a more comfortable transfer. Hand made of fire-proof, water-proof and anti-tear nylon with double seam and outer raised sides to contain patient. Gel sheets are glued in a special "pocket" created to ensure a long duration. **Suitable for all kind of self loading stretcher.**

Dimensions: cm 179x47x13 – weight kg 9,7

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

Length (cm)	179
Width (cm)	47
Thickness (cm)	13
Weight (kg)	9,7
Material	Skay

IMAGES



ACCESSORIES



7018/A-001

Orange skay cushion for mattress
