



## 12080

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher




### FEATURES

The "MBP" platform is a system for Frog and Mercury self-loading stretchers that allows the safe transport of bariatric patients. The stainless steel frame along with the aluminium platform ensure lightness and sturdiness at the same time. The device is completed by belts with reel system on both sides and 4 fixing brackets for anchoring it to the stretcher frame.

**Dimensions: (L x L x H) 1375 x 935 x 130 mm**

**Note:** it is recommended to indicate the stretcher model on which the MPB platform will be installed, so to adjust the fixing clamps.

### MAIN FEATURES

Bariatric	
Length (cm)	137,5
Width (cm)	93,5
Height (cm)	13
Weight (kg)	18
Loading capacity (kg)	350

### IMAGES



---

SEE ALSO

---



**MBS-01 MEBER BARIATRIC  
STRETCHER 7270**  
Self loading stretcher for bariatric  
transport



**12080-001**  
Monting kit for bariatric platform  
on self loading stretcher

---

ACCESSORIES

---



**12080-001**  
Monting kit for bariatric platform  
on self loading stretcher



**7015**  
Folding sky mattress black for  
bariatric MBP platform

ACCESSORY OF



**MBS-01 MEBER BARIATRIC STRETCHER 7270**  
Self loading stretcher for bariatric transport



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite Cinque 7095/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Lite Cinque 7097/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher - high version



**7265 Proof**  
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg



**Mercury 7070 PROOF**  
Certified stainless steel variable height stretcher



**Mercury 7070/4RG PROOF**  
Stainless steel variable height stretcher certified with 4 rotating wheels



**Mercury 7072 PROOF**  
Certified stainless steel variable height stretcher high version



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with variable height and four swivel castors



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**Frog 7210**  
"FROG" yellow self loading stretcher with trendelenburg



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher with trendelenburg certified



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher with trendelenburg certified



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium stretcher with trendelenburg certified