



## 12018N

---

Telescopic foldable phleboclysis stand for self loading stretchers



---

### FEATURES

"P-120" Aluminium telescopic foldable phleboclysis stand for self loading stretchers.

**Dimensions:** min. height cm 45 - max. height cm 128

Height 1st stage: cm 70

Height 2nd stage: cm 98

---

### IMAGES



SEE ALSO



**936**  
Phlebotomy stand for automatic  
loading stretcher Winner and Frog



**936/A**  
Phlebotomy stand for stretchers  
Art. 1510 Proof and Art. 1515 Proof

ACCESSORY OF



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with  
variable height and four swivel  
castors



**Mercury 7070 PROOF**  
Certified stainless steel variable  
height stretcher



**Mercury 7070/4RG PROOF**  
Stainless steel variable height  
stretcher certified with 4 rotating  
wheels



**Mercury 7072 PROOF**  
Certified stainless steel variable  
height stretcher high version



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable  
height stretcher high version with 4  
rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels -  
high version



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version



**MeBer-X 7300**  
"MeBer-X" Multi-level stretcher  
complete with mattress and belts