



1880-001

Mattress for Icarus



FEATURES

Mattress for Icarus made of inner fireproof padding and fireproof, waterproof, anti-tear Skay cover. Double seam.

Dimensions: cm 180x55x5.

*All Pictures shown are for **illustration purpose** only.

MAIN FEATURES

Length (cm)	180
Width (cm)	55
Thickness (cm)	5
Material	Skay

SEE ALSO



1880-002

Mattress pillow to be used with
Icarus stretcher

ACCESSORIES



Icarus 13010
Multi level stretcher with side rails



1880-002
Mattress pillow to be used with
Icarus stretcher

ACCESSORY OF



Icarus 13010
Multi level stretcher with side rails
