




19140

"MTX-2" double space stretcher support with lateral movement and ramp of mm. 350 - EN1789 certified




FEATURES

Double space stretcher support with lateral movement and ramp mm 350 in the underneath space can be stored both a scoop stretcher and a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956.
Gas spring included.



Lateral movement: cm 30

MAIN FEATURES

Side sliding	
Width (cm)	63
Length (cm)	229
Weight (kg)	76
Loading capacity (kg)	400
Plate height	25
Total height	28

IMAGES



ACCESSORIES



19110
Stainless steel plate mm 1917 x 572
for stretcher supports



787L
Gas spring for ramp



19123
Ribaltina 130 mm con attacchi per
supporti e rullo



19115
Ribaltina da 180 mm con attacchi
per supporti e rullo

STANDARD ACCESSORIES



821G
Gas spring for ramp