



## 1080

Rescue mask



### FEATURES

Indispensable for all rescue operators, this mask avoid any contacts between mouths during emergency reanimations. Made of transparent material, it is composed by a soft cushion and a valve with non-return membrane. It is all contained in a handy, plastic, shockproof case. Color of case may change.

### MAIN FEATURES

Size	Adults
Material	PVC

### SEE ALSO



**1787**  
Keyring with mouth-to-mouth  
mask

---

ACCESSORY OF

---



**1366**

Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---