



2401


Silicone reanimation mask size 1



FEATURES

Reanimation and anesthesia mask. Anatomical border in silicone guarantees perfect fitting to face while the shell, made in transparent material, allows the patient to be monitored during the application. Standard connector. **Can be autoclave sterilised, available in 6 sizes.**

MAIN FEATURES

Size	Children
Pediatric	
Material	Silicone

SEE ALSO



2402

Silicone reanimation mask size 2



2401

Silicone reanimation mask size 1



2404

Silicone reanimation mask size 4



2405

Silicone reanimation mask size 5



2400

Silicone reanimation mask size 0



2403

Silicone reanimation mask size 3

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE