



417

“NEOFIX” Neoprene splints set with case



FEATURES

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanification. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	

---

SEE ALSO

---



**Halley 886**  
Vacuum splint set with pump and carrying case



**Neofix 412**  
Neoprene splint for wrist



**Neofix 410**  
Neoprene splint for forearm



**Neofix 415**  
Neoprene splint for ankle/elbow



**Neofix 407**  
Neoprene splint for arm



**Neofix 405**  
Neoprene splint for leg