



## 616/A

Certified vertical fixing system for scoop stretchers



### FEATURES

Aluminium supports with safe and rapid hook to fix the scoop stretcher in vertical position to the wall or inside the ambulance. With our item 616/A it is possible to fix the scoop stretcher folded.



Device certified according to the European harmonized safety standards UNI EN 1789.

### MAIN FEATURES

Length (cm)	44
Width (cm)	2,5
Height (cm)	9,5
Weight (kg)	0,55


### SEE ALSO



**616**

Certified vertical fastening system  
for scoop stretchers

ACCESSORIES



**Maxima 630/V**  
Green painted anti trauma scoop  
stretcher



**Maxima 630/G**  
Yellow painted anti trauma scoop  
stretcher



**Maxima 630/B**  
Blue painted anti trauma scoop  
stretcher



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher