



# 630/G Maxima

---

Yellow painted anti trauma scoop stretcher



---

## FEATURES




Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

## MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	
Material	Painted aluminium
Trauma	
Alpine	

## SEE ALSO



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher



**Maxima 630/N**  
Black painted anti trauma scoop  
stretcher



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher



**Maxima 630/B**  
Blue painted anti trauma scoop  
stretcher



**Maxima 630/V**  
Green painted anti trauma scoop  
stretcher



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc

## ACCESSORIES




**616**  
Certified vertical fastening system  
for scoop stretchers



**614**  
Certified horizontal fastening  
system for scoop stretchers



**630/CS**  
Carrying case for scoop stretchers



**616/A**  
Certified vertical fixing system for  
scoop stretchers



**603/MEB**  
Nylon bag containing 3 belts Art.  
607/Ameb