




630/V Maxima

Green painted anti trauma scoop stretcher



FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
Type	
10G Certification	<input checked="" type="checkbox"/>
Material	Painted aluminium
Trauma	<input checked="" type="checkbox"/>
Alpine	<input checked="" type="checkbox"/>

SEE ALSO



Maxima 630/N
Black painted anti trauma scoop
stretcher



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630/A
Orange painted anti trauma scoop
stretcher



Maxima 630/G
Yellow painted anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 630
Grey anodized anti trauma scoop
stretcher

ACCESSORIES



630/CS

Carrying case for scoop stretchers



616/A

Certified vertical fixing system for
scoop stretchers



616

Certified vertical fastening system
for scoop stretchers



614

Certified horizontal fastening
system for scoop stretchers



603/MEB

Nylon bag containing 3 belts Art.
607/Arneb