



## 643 Folding stretcher


Collapsible rescue stretcher foldable in 4 pieces



### FEATURES

Built in tubular anodized aluminium, canvas in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Ideal for storage inside of any emergency vehicle. Straps are not included. (see "Restraint systems" Art. 604/MEB).

### MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Min height (cm)	12
Width (cm)	49
Max length (cm)	195
Folded length (cm)	97
Epaisseur pliée (cm)	17
Weight (kg)	5
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

### IMAGES



---

## SEE ALSO

---



**Folding stretcher 646**  
Collapsible rescue stretcher  
foldable in 2 pieces

---

---

## ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon

---