



866 Sirius

Rigid cervical collar in 1 piece size STANDARD



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	

SEE ALSO



Antares 854
Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Antares 856
Cervical collar set 6 pieces with carrying case



Sirius 874
Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Sirius 876
Cervical collar set 6 pieces with carrying case



Sirius 873
Cervical collar set 3 pieces with carrying case Art. 864 - 866 - 868

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE
