



9022 Duck

Yellow pediatric spinal board



FEATURES

Designed for immobilization and transport of children with height up to 120 cm. This spinal board has been produced with plastic material at high strenght and studied paying particular attention to the problems connected to daily interventions. It has 4 holes for the quick and total fixing of the head immobilizer and 2 cavities where the board lays on the floor, when base is blocked in traditional way, that allow to avoid damages to the rip-off straps during the usage of the spinal board or its accomodation in the ambulance.

It has 10 handles for the transport.

It allows preliminary disgnostic under X-ray and it's available in 3 different colours: yellow, orange and blue.
It isn't supplied with belts.

MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	125
Width (cm)	41
Thickness (cm)	5
Weight (kg)	3,7
Loading capacity (kg)	100
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	ABS

SEE ALSO



Duck 9024
Orange pediatric spinal board



Duck 9026
Blue pediatric spinal board

ACCESSORIES



Kids 627
Padded head immobilizer –
pediatric



Spider Baby 692
Immobilization belts for pediatric
spinal board



16150-002
Storage bag for Ergon stretcher