



## 938

Oximeb oxygen cylinder holder



### FEATURES

Cylinder holder produced with strong, antitear nylon and inside padding to prevent damages from knocks. It is supplied with a handle for the transport and a side pocket for tube and mask.

Can be used with cylinders of 2 liters and, thanks to a handy fastening, it can be mounted on the side protection of the automatic loading stretchers and on the back rest of the sedan chairs.

### MAIN FEATURES

Weight (kg)	0,3
Length (cm)	31
Width (cm)	19
Depth (cm.)	13
Material	Nylon

### SEE ALSO



12054

Oximeb Maxi - Oxygen cylinder holder

---

ACCESSORY OF

---



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg