



939

Instrument support for stretchers



FEATURES

Instrument support to be fixed on Me.Ber. stretcher frames' to transport medical equipments, defibrillators and monitors. It is supplied with belts to fix devices and is foldable to be stored inside ambulances.

Dimensions: open cm 34x56x37h, closed height cm 8

MAIN FEATURES

Width (cm)	56
Depth (cm.)	34
Height (cm)	37
Epaisseur pliée (cm)	8
Weight (kg)	5,6
Loading capacity (kg)	30

IMAGES



SEE ALSO



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels

ACCESSORY OF



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg