



955 PROOF

EN 1789 fixing system for Proof roll-in stretchers



FEATURES

Designed to guarantee maximum stability of the stretcher on board the rescue vehicle, they are made entirely of stainless steel for longer life and corrosion resistance. Their easy-to-install structure ensures that the stretcher or trolley, and patient are securely restrained in all directions. Art. 955 Proof is suitable for fixing all our 10G. roll-in stretchers whilst art. 956 Proof is suitable for fixing "Titanus" trolleys.



Device certified according to the European harmonized standards UNI EN 1789.

Height	Width	Lenght	Weight
13 cm	48 cm	177 cm	7,7 kg

MAIN FEATURES

Height (cm)	13
Length (cm)	177
Weight (kg)	7,7

SEE ALSO



956 PROOF

Fixing system for Titanus automatic loading trolley

ACCESSORIES



955-001

Front reverse-plate for 955-956 Proof certified fixing system



955-002

Rear reverse-plate for 955-956 Proof certified fixing system

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270
Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with variable height and four swivel castors



Mercury Lite 7082/4RG Proof
Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury 7070 PROOF
Certified stainless steel variable height stretcher



Mercury 7070/4RG PROOF
Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072 PROOF
Certified stainless steel variable height stretcher high version



Mercury 7072/4RG PROOF
Certified stainless steel variable height stretcher high version with 4 rotating wheels



Frog 7240 PROOF
"FROG PLUS" self loading stretcher with trendelenburg certified



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium stretcher with trendelenburg certified



Frog 7210 PROOF
"FROG" yellow self loading stretcher with trendelenburg certified



Mercury Lite Cinque 7095/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite 7082/4RG Proof
Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury EL 7094 PROOF
Certified self-loading stretcher - high version



7265 Proof
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg