



## 425

Arcos long rigid splint for lower limbs



### FEATURES


New rigid splints, made with fire-proof and water proof PVC, allow a correct immobilization of the limbs. They are X-ray translucent. Equipped with 4 Velcro® closures with adjustable lenght allowing operator to apply them easily and safely. Shaped to strongly and easily fix to ankle and upper part of the leg, they are supplied with foot support integrated.

Available in two sizes:


Adult arm/ pediatric leg: lenght cm. 52 - base cm 9

Adult leg: lenght cm 74 - base cm 13

### MAIN FEATURES

Trauma	
Width (cm)	27
Length (cm)	74
Weight (kg)	1,1
Usage	Immobilization
Material	PVC

### IMAGES



---

## SEE ALSO

---



**430**

Arcos rigid splint set with carrying bag



**427**

Arcos short rigid splint for upper limbs

---