



9124/B


Standard spinal immobilizer, blue colour



FEATURES

Essential device used to safely immobilize and extricate a patient. Made of vinyl-coated nylon. Complete with a padded pillow, two head straps and carrying bag.

MAIN FEATURES

Trauma	
Usage	Immobilization
Width (cm)	85
Height (cm)	83
Belts' length (cm)	120
Max circumference	180
Weight (kg)	2,9
Loading capacity (kg)	230
Material	Nylon

SEE ALSO



Med 899
Spinal immobilizer



Med Ped 9120
Pediatric spinal immobilizer