



## 678 Extra Lite Basic

Extra Lite Basic, foldable carrying chair



### FEATURES

EXTRA LITE basic is the new Me.Ber. stair chair entirely made of lightweight aluminium alloy frame. The backrest and back carrying handles inclination, along with the unique frontal arms configuration, enable a perfect weights distribution, relieving operators' efforts. While manoeuvring on smooth and plane surfaces, the 200 mm rear and 125 mm front wheels gives exceptional comfort. The structure along with an ergonomic seat and backrest, complete the characteristics of this new rescue device.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

### MAIN FEATURES

Type	
Foldable	
Track assisted system	
Height (cm)	101
Width (cm)	55
Length (cm)	76
Max length (cm)	139
Epaisseur pliée (cm)	29
Weight (kg)	13,6
Loading capacity (kg)	230
Folded height (cm)	101
Cardio	
10G Certification	
Material	Aluminium alloy

---

## SEE ALSO

---



**679**  
Narrow foldable sedan chair



**Extra Lite Standard 676**  
Foldable carrying chair with 4 wheels, footrest and arms

---

## ACCESSORIES

---



**670-001**  
Orange/black head immobilizer for Extra stair chair



**960**  
"Extra fixing system" - fixing system for chair Extra



**566**  
Black bag with 3 belts in 2 pieces made of black ribbon with orange MeBer logo for Extra chairs



**975**  
"Ergolift fixing system" - Fixing system for chair Extra Ergolift



**670-004**  
Castor with brake diam. 200 mm for Extra chair



**670-002**  
Transpire, dust proof black nylon cover for Extra chairs

---

STANDARD ACCESSORIES

---



566

Black bag with 3 belts in 2 pieces  
made of black ribbon with orange  
MeBer logo for Extra chairs

---