



TECHNICAL DATASHEET

Item code	7070/4RG PROOF	
Name	MERCURY - Stainless steel certified variable height self-loading stretcher with 4 rotating wheels	
Dimensions [mm]	Length	1970
	Length with folded trolley	1650
	Width	570
	Max height	950
	Intermediate height	720
	Folded height	310
	Loading height	640 ÷ 660
Used materials	AISI 304, ABS, HDPE	
Wheels	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm	
Weight [kg]	44	
Max load [kg]	250	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
Conformity	MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it