



## TECHNICAL DATASHEET

<b>Item code</b>	<b>7070 PROOF</b>	
<b>Name</b>	<b>MERCURY - Stainless steel self loading stretcher with variable height</b>	
<b>Dimensions [mm]</b>	Length	1970
	Length with folded trolley	1650
	Width	570
	Max height	950
	Intermediate height	720
	Folded height	310
	Loading height	640 ÷ 660
<b>Used materials</b>	AISI 304, ABS, HDPE	
<b>Wheels</b>	n°4 Ø200 mm (2 fixed+2 swivel) n°2 Ø150 mm n°2 Ø60 mm	
<b>Weight [kg]</b>	42	
<b>Max load [kg]</b>	250	
<b>Declared use</b>	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device EN 1865-1 - EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)