



TECHNICAL DATASHEET

| | | |
|------------------------|---|-----------|
| Item code | 7072 PROOF | |
| Name | MERCURY - Stainless steel certified high version self loading stretcher with variable height | |
| Dimensions [mm] | Length | 1970 |
| | Length with folded trolley | 1650 |
| | Width | 570 |
| | Max height | 990 |
| | Intermediate height | 750 |
| | Folded height | 310 |
| | Loading height | 690 ÷ 710 |
| Used materials | AISI 304, ABS, HDPE | |
| Wheels | n°4 Ø200 mm (2 fixed+2 rotating) n°2 Ø150 mm n°2 Ø60 mm | |
| Weight [kg] | 43 | |
| Max load [kg] | 250 | |
| Declared use | Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled. | |
| Conformity | MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified | |



2204

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it