



TECHNICAL DATA SHEET

Item code	7082/4RG PROOF	
Name	MERCURY LITE – Aluminium alloy self loading stretcher with variable height with 4 swivel wheels high version	
Dimensions [mm]	Length	1970
	Length with folded trolley	1650
	Width	570
	Max height	990
	Intermediate height	750
	Folded height	310
	Loading height	690 ÷ 710
Used materials	Aluminium alloy, HDPE	
Wheels	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm	
Weight [kg]	40,5	
Max load [kg]	250	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled.	
Conformity	MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it