

# MEBER+

IN CASE OF EMERGENCY

**www.merber.it**

Company press service realized by Me.Ber. s.r.l.  
For information please, contact:  
contact@meber.it or go to the link: www.merber.it

**Medical Journal**

## EXTRA<sup>ERGO</sup>LIFT

POWERED STAIR CHAIR

Me.Ber. presents:  
**EXTRA ERGOLIFT**  
the revolutionary  
powered stair  
chair

ORDINARY.  
EXTRAORDINARY.  
FOCUS ON:  
**EXTRA ERGOLIFT**

**3**

WHY CHOOSE  
**EXTRA ERGOLIFT?**

**4**

NEWS

**8**

18 ramps of 20 steps  
with a load of 130 kg



**TRAUMA**

Aluminum alloy

Discover it! >>>

**Upward to an higher level**




Use this QR Code to visit Extra Ergolift page on our web site.

# UPWARD TO AN HIGHER ...LEVEL

Meber presents its first revolutionary motorized Extra Ergolift chair to ascend and descend stairs. It derives from the successful Extra series, maintaining the main ergonomic and structural features: innovative design, wide use of Hi-Tech materials, extreme easiness of use, lightness and reliability.



*Extra Ergolift*



*Extra Ergolift on fixyng system*

The main fields of application for a powered chair are typically the ordinary and extraordinary services, for example the transport services to pick up a patient or take him back

home or to a private hospital, where an elevator may not be present. Day by day, while working on the Extra Ergolift project, we realized that its unique features enable it also to

be used in critical situations or as a support patients handling device in large structures (airports, stations, hotels, etc).





# ORDINARY. EXTRAORDINARY.

We are proud to present **Extra Ergolift** as **ordinary** and **extraordinary** device. Rescue crews consisting of women or weak people are able to handle even heavy patients without any effort.

**Extra Ergolift** allows also a safe transport of bariatric patients up to 200 kg weight in motorized transport configuration and 230 kg in manual, both upward and downward stairs.



### User friendly

➔ Another proof of the **usefulness** and of the **strenght** of this product is the possibility to retrieve from basements no-traumatic patients fastly and safely.

With **Extra Ergolift** a crew, composed by two people only, can evacuate a patient without **any difficulty**.

➔ **Extra Ergolift** is the result of an intense course of research, innovation and improvement, three important values for Me.Ber.

Extra Ergolift is the connection of the most innovative structural and ergonomic expedients. The best of innovation shaped the most advanced and safest chair in Me.Ber.



# WHY CHOOSE EXTRA ERGOLIFT?

## THE ANSWER IN 6 POINTS:

Unique device for everyone! For slender and strong people, for old and young, everyone can handle the powered Extra Ergolift chair with minimum effort. It doesn't matter the type of stair: linear, external or spiral staircase. Extra Ergolift has unique features which allow also the transport of bariatric people in every circumstance.

**WIRELESS CONSOLE**

**LIGHT ALUMINIUM  
ALLOY**

**SIMPLE AND  
COMPACT ENGINE**

**HIGH COMFORT  
LEVEL**

**LONG LIFE  
BATTERIES**


**STABILITY ON  
STAIRS**





# STRENGTH POINTS

## 1. LONG AND HIGH-GRIP SLIDING TRACKS




The **Extra Ergolift** sliding tracks have a longer length allowing to constantly lay on three steps for better stability and comfort.

## 2. UNMATCHED STABILITY AND SAFETY

Everyone who choose **Extra Ergolift** has a further security: to guarantee the maximum level of safety as required by law. Active and passive safeties ensure a safer work even in critical situations.

**Extra Ergolift** is designed to prevent injuries during the service and to guarantee maximum efficiency even for weak operators. Our chair allows maximum efficiency in all conditions, ensuring a safe and comfortable transport for patients.



### 3. HIGH AUTONOMY



The engine is simple and compact while the battery has a long life (Milwaukee M28 RED LITHIUM 5,0 Ampere per Extra Ergolift chair.)

### 4. IP65 ISOLATION CODE

The IP Code, International Protection Marking, classifies and rates the degree of protection provided against intrusion (body parts such as hands and fingers), dust, accidental contact, and water by mechanical casings and electrical enclosures.

Tested with IP 65 index protection from water and dust, Extra Ergolift can be used under the rain or with adverse weather conditions.



The prefix IP is followed by two numbers:

- The first digit indicates the level of protection that the enclosure provides against access to hazardous parts (e.g., electrical conductors, moving parts) and the ingress of solid foreign objects.
- The second digit indicates the level of protection that the enclosure provides against harmful ingress of water.

**IP6X= Totally protected against dust and every little solid objects; IPX5= Protected against water gush**



## NOT TO MISS



**For the transport  
of your chair in  
ambulance, choose  
the appropriate  
fixing system:  
Art. 975 Extra  
Ergolift fixing  
system.**

## MeBer thinks of You. Always.

To manage the organization of the crews in order to find right combination between available rescue operators according to their muscular strength, is very complicated.

Likewise to manage cases of injuries of rescuers during their duties, maybe for wrong movements or not correct operations.

The optimal management of the crews is an important topic in order to guarantee the maximum safety, not only for patients, but also for rescuers.

Extra Ergolift has been designed to avoid injuries during the service and to guarantee the maximum operativity, also to crews with "weak people". Everyone can handle heavy patients without any effort on the stairs.

## SUPPLIED WITH:



**Black bag with three belts  
in two pieces**



**Battery charger + Battery  
MILWAUKEE 5,0 Ampere**

## ACCESSORIES



**Transpire cover for Extra  
Ergolift**

# STAY TUNED! WE ARE ALSO ONLINE



## Short news

Me.Ber.: thirty years in the emergency world

2017 is an important year for Me.Ber.: 30 years of career in the world of emergency. Our aim has always been to give quality tools to rescuers and medical emergency workers. We dedicate to them our new product catalogue cover, because everyone brought us ideas, support and above all, trust. We would like to thank one by one and also share with you our journey towards innovation, because Me.Ber. is also you! So happy birthday!

Join us on our social networks and subscribe to our newsletter

Fairs, congress, events, promotion, news and much more. We will update you about every company news and not only. Go to our website: [www.meber.it](http://www.meber.it).

Me.Ber. Promo/used Catalogue 2017

Do you already know our promo/used catalogue? Click [here](#) to discover all the discounted products or call: +39/ 0521.648770.

