

MEBER



MEBER-X MULTILEVEL STRETCHER



Art. 7300 Accessories not included



MEBER



Adjustable backrest:
0÷80°

Weight: 45,9 kg
MAX Load: 320 kg

Trendelenburg: 15°
Fowler: 35°

Length: 1660 ÷ 1970
Width: 600

4 rotating
wheels

Rear wheels
with brake

Materials:
Aluminium
alloy
Polyamide



MEBER-X MULTI LEVEL STRETCHER

HEIGHT		LOADING HEIGHT	
POSITION 1	385	POSITION 1	0
POSITION 2	540	POSITION 2	180
POSITION 3	640	POSITION 3	290
POSITION 4	730	POSITION 4	390
POSITION 5	810	POSITION 5	485
POSITION 6	885	POSITION 6	580
POSITION 7	955	POSITION 7	670
POSITION 8	1020	POSITION 8	760
POSITION 9	1080	POSITION 9	855