



894/SNOW Snake Snow

Independent chambers vacuum mattress with anti-slip system for snow



FEATURES

Vacuum mattress with 14 independent chambers for an excellent and fast immobilisation. **Complete with anti-slip system for rescue on snow and ice.** Extremely light and sturdy, 8 ergonomics handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Usable between -30°C and +70°C. An unmatched vacuum durability without any deflection, **the UNI EN1865 certificate**, according to the European harmonized standards, **complete features of this device.**

MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	210
Width (cm)	100
Weight (kg)	8,9
Loading capacity (kg)	150
Trauma	<input checked="" type="checkbox"/>
Alpine	<input checked="" type="checkbox"/>
Material	Polyurethane

IMAGES



SEE ALSO



Chesterfield 9114
Vacuum mattress with
independent chambers with
central welding

ACCESSORIES



888
Aluminium aspiration pump



896/DF
Abs aspiration pump
