



894/SNOW Snake Snow



Independent chambers vacuum mattress with anti-slip system for snow



FEATURES

Vacuum mattress with 14 independent chambers for an excellent and fast immobilisation. **Complete with anti-slip system for rescue on snow and ice.** Extremely light and sturdy, 8 ergonomics handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Usable between -30°C and +70°C. An unmatched vacuum durability without any deflection, **the UNI EN1865 certificate**, according to the European harmonized standards, **complete features of this device.**

MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	210
Width (cm)	100
Weight (kg)	8,9
Trauma	
Alpine	
Material	Polyurethane
Loading capacity (kg)	150

IMAGES



SEE ALSO



Chesterfield 9114
Vacuum mattress with
independent chambers with
central welding

ACCESSORIES



888
Aluminium aspiration pump



896/DF
Abs aspiration pump
