



# 894/SNOW Snake Snow

Independent chambers vacuum mattress with anti-slip system for snow



#### **FEATURES**

Vacuum mattress with 14 indipendent chambers for an excellent and fast immobilisation. **Complete with anti-slip system for rescue on snow and Ice.** Extremely light and sturdy, 8 ergonomics handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Usable between -30°C and +70°C. An unmatchable vacuum durability without any deflection, **the UNI EN1865 certificate**, according to the European harmonized standards, **complete features of this device**.

#### MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	210
Width (cm)	100
Weight (kg)	8,9
Loading capacity (kg)	150
Trauma	<b>₩</b>
Alpine	<b>₩</b>
Material	Polyurethane

#### **IMAGES**





## **SEE ALSO**



# Chesterfield 9114 Vacuum mattress with independent chambers with central welding

## **ACCESSORIES**



**888**Aluminium aspiration pump



**896/DF**Abs aspiration pump