



## 1274

---

Silicone reanimation bag for adults



---

### FEATURES

Silicone reanimation bag. Transparent, autoclavable, complete with safety valve with adjustment at 60 mbar for adults. Equipped with connection for oxygen bottle, reservoir and masks.

---

### MAIN FEATURES

Size

Adults

---

### SEE ALSO



**1276**

Silicone reanimation bag for children



**1278**

Silicone reanimation bag for newborns

---

ACCESSORY OF

---



**1366**  
Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---