



11282 Clineb Basic

Professional piston compressor nebulizer for aerosoltherapy



FEATURES

CLINEB BASIC is a piston-type compressor nebulizer system working at 230V/50Hz (other voltages available upon request). High performance with any type of drug, ideal for intensive hospital and clinic use. Manufactured with high thermal and electric insulation plastic chassis in compliance with the latest European Safety regulations. The oil-free piston compressor has long durability and is equipped with the highly efficient HI-FLO jet nebulizer to guarantee quick and accurate drug delivery. The device is designed for easy transport and handling and is recommended for atomising antibiotics and bronchodilator drugs. The medical device is designed for continuous use.

STANDARD EQUIPMENT: • HI-FLO Jet Nebulizer • Adult Mask • Pediatric Mask • Mouthpiece • Air Tube • Air Filter