



13012 Icarus

Multi level stretcher with steel tank




FEATURES

Multi level stretcher Icarus is made with stainless steel, has 4 wheels 2 of which with brakes, side rails and foot rest. Four gas pistons allow a perfect adjustment to all necessary heights for a practical and professional work.

Available in following versions with:

- ABS lying part for a perfect cleaning and disinfection (art. 13010)
- stainless steel tank (art. 13012)

MAIN FEATURES

Multilevel	
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Min height (cm)	51
Max height (cm)	96
Width (cm)	62
Length (cm)	198
Weight (kg)	39
Loading capacity (kg)	170
Size	Adults
Usage	Transport
Material	Stainless steel

SEE ALSO



Icarus T3010

Multi level stretcher with side rails

ACCESSORIES



607/B

Blue belt in 2 pieces with plastic buckle
