



7260 proof Frog Lite 7260 PROOF

"FROG LITE" self loading aluminium stretcher with trendelenburg certified



FEATURES

Made of painted aluminium alloy, extremely light, it is also provided with trendelenburg. The 4 Ø 200 mm polyurethane wheels (2 fixed with brake and 2 rotating with self-centring system) grant high comfort on uneven surfaces. Easy to clean and disinfect.

- Rigid platform made of painted aluminium alloy
- Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- Folding sides
- Vertical legs protected by plastic buffets
- Rotating wheels are self centring
- Collapsible frontal trolley


In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

10G Certification	
Height (cm)	91,5
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	31
Width (cm)	56
Length (cm)	196
Weight (kg)	36,9
Loading capacity (kg)	250
Filter feature	
Material	Painted aluminium HDPE (High-density polyethylene)

IMAGES



SEE ALSO



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified



Frog 7210 PROOF
"FROG" yellow self loading stretcher
with trendelenburg certified

ACCESSORIES



955 PROOF
EN 1789 fixing system for Proof roll-
in stretchers



939
Instrument support for stretchers



936
Phleboclysis stand for automatic
loading stretcher Winner and Frog



12025
Black multi usage object holder for
automatic loading stretcher



12106
Mattress for AIR-COM stretchers -
Anti-bedsore mattress with
alternate segments



7015
Folding sky mattress black for
bariatric MBP platform



12154
"BEAR PODS" Bariatric
containment system for self-
loading stretchers



707
Red 6-handle fireproof cover with
Me.Ber tape for self-loading
stretcher



938
Oximeb oxygen cylinder holder



12054
Oximeb Maxi - Oxygen cylinder
holder



12128
Adjustable system for headrest of
Meber-X and Mercury stretchers



12080
"MPB" wide bariatric platform for
self-loading stretchers Frog
Mercury and Gecko complete of
belts with reel and fixing for
stretcher



12080-001
Monting kit for bariatric platform
on self loading stretcher



12070/BG
"COVER OX 2-3 BG"- professional
cylinder holder for 2 and 3 lt oxygen
bottles



12070/BR
COVER OX2-3 BR professional
cylinder holder with belt and
orange MeBer logo for 2 and 3 lt
oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR

"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles



12084

Bariatric wings set for self-loading
stretchers complete with
mattresses and belts with reel

STANDARD ACCESSORIES



7013/N

3 pcs. Skay fire-proof electrowelded
mattress



612/N-2/MEB

Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



608/N-MEB

2 pieces black belt with metalic
hook and MeBer band (L. 155 cm)