



7260 proof Frog Lite 7260 PROOF

"FROG LITE" self loading aluminium stretcher with trendelenburg certified



FEATURES

Made of painted aluminium alloy, extremely light, it is also provided with trendelenburg. The 4 Ø 200 mm polyurethane wheels (2 fixed with brake and 2 rotating with self-centring system) grant high comfort on uneven surfaces. Easy to clean and disinfect.

- Rigid platform made of painted aluminium alloy
- Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- Folding sides
- Vertical legs protected by plastic buffets
- Rotating wheels are self centring
- Collapsable frontal trolley


In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

10G Certification	
Height (cm)	91,5
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	31
Width (cm)	56
Length (cm)	196
Weight (kg)	36,9
Loading capacity (kg)	250
Filter feature	
Material	Painted aluminium HDPE (High-density polyethylene)

IMAGES



SEE ALSO



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified



Frog 7210 PROOF
"FROG" yellow self loading stretcher
with trendelenburg certified


ACCESSORIES



939
Instrument support for stretchers



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



12048
Certified interface for incubator fixation for self loading trolleys



12025
Black multi usage object holder for automatic loading stretcher



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



12106
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



7015
Folding sky mattress black for bariatric MBP platform



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



707
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



938
Oximeb oxygen cylinder holder



12054
Oximeb Maxi - Oxygen cylinder holder



12128
Adjustable system for headrest of Meber-X and Mercury stretchers



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001
Monting kit for bariatric platform on self loading stretcher



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

STANDARD ACCESSORIES



7013/A
3 pcs. Skay fire-proof electrowelded
mattress



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)
