



## 7260 proof Frog Lite 7260 PROOF

---

"FROG LITE" self loading aluminium stretcher with trendelenburg certified



---

### FEATURES

Made of painted aluminium alloy, extremely light, it is also provided with trendelenburg. The 4 Ø 200 mm polyurethane wheels (2 fixed with brake and 2 rotating with self-centring system) grant high comfort on uneven surfaces. Easy to clean and disinfect.

- Rigid platform made of painted aluminium alloy
- Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- Folding sides
- Vertical legs protected by plastic buffets
- Rotating wheels are self centring
- Collapsable frontal trolley

In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB).




Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

---

## MAIN FEATURES

---

|                         |   |
|-------------------------|---|
| 10G Certification       |  |
| Height (cm)             | 91,5  |
| Max loading height (cm) | 66  |
| Min loading height (cm) | 60  |
| Folded height (cm)      | 31  |
| Width (cm)              | 56  |
| Length (cm)             | 196   |
| Weight (kg)             | 36,9  |
| Loading capacity (kg)   | 250   |
| Filter feature          |   |
| Material                | Painted aluminium   |

---

## IMAGES

---



---

## SEE ALSO

---



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher  
with trendelenburg certified

---

## ACCESSORIES



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**939**  
Instrument support for stretchers



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**12025**  
Black multi usage object holder for automatic loading stretcher



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



**7015**  
Folding sky mattress black for bariatric MBP platform



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**707**  
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



**938**  
Oximeb oxygen cylinder holder



**12128**  
Adjustable system for headrest of Meber-X and Mercury stretchers



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**  
"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**12084**

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

---

---

## STANDARD ACCESSORIES

---



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)

---