



## 7260 proof Frog Lite 7260 PROOF

---

“FROG LITE” self loading aluminium stretcher with trendelenburg certified



---

### FEATURES

Made of painted aluminium alloy, extremely light, it is also provided with trendelenburg. The 4 Ø 200 mm polyurethane wheels (2 fixed with brake and 2 rotating with self-centring system) grant high comfort on uneven surfaces. Easy to clean and disinfect.

- Rigid platform made of painted aluminium alloy
- Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- Folding sides
- Vertical legs protected by plastic buffets
- Rotating wheels are self centring
- Collapsable frontal trolley

In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB).




Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

---

## MAIN FEATURES

---

10G Certification	
Height (cm)	91,5
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	31
Width (cm)	56
Length (cm)	196
Weight (kg)	36,9
Loading capacity (kg)	250
Filter feature	
Material	Painted aluminium HDPE (High-density polyethylene)

---

## IMAGES

---



---

## SEE ALSO

---



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher  
with trendelenburg certified

---

## ACCESSORIES



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**939**  
Instrument support for stretchers



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**12025**  
Black multi usage object holder for automatic loading stretcher



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



**7015**  
Folding sky mattress black for bariatric MBP platform



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**707**  
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



**938**  
Oximeb oxygen cylinder holder



**12054**  
Oximeb Maxi - Oxygen cylinder holder



**12128**  
Adjustable system for headrest of Meber-X and Mercury stretchers



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**12084**

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

---

## STANDARD ACCESSORIES

---



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)

---