



## 644/S Rigid stretcher

Stiff fixed stretcher for storage with backrest



### FEATURES

Fixed rigid stretcher for storage, featuring a 3 positions adjustable backrest. Supplied with 2 quick-release straps for patient immobilization art. 604/MEB.

### MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Height (cm)	17
Width (cm)	48
Max length (cm)	186
Weight (kg)	5,9
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

### SEE ALSO



**Fixed stretcher 644**  
Stiff fixed stretcher for storage

---

## STANDARD ACCESSORIES

---



**604/MEB**

Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon

---