



11010/FA Cormos

Portable aspirator for secretions with autoclavable canister



FEATURES

Portable aspirator for secretions CORMOS with autoclavable liner, green MeBer bag, 12V battery charger, suction rigid probe with hose and rechargeable bracket 10G. Indicated for emergency services. CORMOS works directly with the 12V. cc. of the vehicle or with an inner accumulator with more than 45 min. of autonomy. CORMOS is silent and powerful thanks to its engine of 24 litres/min which enables a constant aspiration, with control manometer, positioned on the front of the aspirator, superior to 850 mbar. Easy to carry thanks to its particular and compact dimension, the anatomic handle and its limited weight of 4,6 Kg. New case design, made of ABS/UL high resistance anti-shock with dimension of cm 34x11x24. It is supplied with 1000 ml autoclavable canister with safety valve and antibacterial filter to fix directly onto the cap. The "CORMOS" aspirator for secretions is manufactured according to the ISO 10079-1:2009 product normative, IIb class, EN 1789 and is certified CE 0123, according to the directive 93/42/CEE for medical devices and european directive 104 (automotive).

MAIN FEATURES

Width (cm)	34
Height (cm)	24
Depth (cm.)	11
Weight (kg)	4,6

SEE ALSO



Cormos 11010/LN
Portable aspirator for secretions
with disposable liner

ACCESSORIES



11025

Serres disposable liners 1000 ml (x
12 pz./pcs.)



11024

Serres disposable liners 1000 ml (x 6
pz./pcs.)



11034

Battery charger for Cormos 2 poles



11023

10g certified wall bracket with
integrated charging system for
CORMOS and CORMOS PLUS

STANDARD ACCESSORIES



11014

Suction rigid probe and hose for
"CORMOS" and "CORMOS PLUS"
aspirators



11016

Green carrying bag for suction unit