



# **11010/FA** Cormos

Portable aspirator for secretions with autoclavable canister



#### **FEATURES**

Portable aspirator for secretions CORMOS with autoclavable liner, green MeBer bag, suction rigid probe with hose and rechargeable bracket 10G. Indicated for emergency services. CORMOS works directly with the 12V. cc. of the vehicle or with an inner accumulator with more than 45 min. of autonomy. CORMOS is silent and powerful thanks to its engine of 25 litres/min which enables a constant aspiration, with control manometer, positioned on the front of the aspirator, superior to 850 mbar. Easy to carry thanks to its particular and compact dimension, the anatomic handle and its limited weight of 4,6 Kg. New case design, made of ABS/UL high resistance anti-shock with dimension of cm 34x11x24. It is supplied with 1000 ml autoclavable canister with safety valve and antibacterial filter to fix directly onto the cap. The "CORMOS" aspirator for secretions is manufactured according to the ISO 10079-1:2009 product normative, Ilb class, EN 1789 and is certified CE 0123, according to the directive 93/42/CEE for medical devices and european directive 104 (automotive).

### MAIN FEATURES

Width (cm)	34
Height (cm)	24
Depth (cm.)	11
Weight (kg)	4,6

## **SEE ALSO**



Cormos 11010/LN
Portable aspirator for secretions
with disposable liner

#### **ACCESSORIES**



11025 Serres disposable liners 1000 ml (x 12 pz/pcs.)



11024 Serres disposable liners 1000 ml (x 6 pz/pcs.)



**11026** Serres disposable liners 1000 ml (x 36 pz/pcs.)



11162 Disposable 1-litre canister



11012 Autoclavable canister



11013 Antibacterical filter



11022-003 2 poles cable for aspirator CORMOS



11034 Battery charger for Cormos 2 poles

## STANDARD ACCESSORIES



11023
10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS



11028 Autoclavable canister



11014
Suction rigid probe and hose for
"CORMOS" and "CORMOS PLUS"
aspirators



**11016**Green carrying bag for suction unit