



## 11010/LN Cormos

Portable aspirator for secretions with disposable liner



### FEATURES

Medical aspirator for secretions CORMOS, compact and powerful, indicated for emergency services. CORMOS works directly with the 12V. cc. of the vehicle or with an inner accumulator with more than 45 min. of autonomy. CORMOS is silent and powerful thanks to its engine of 24 litres/min which enables a constant aspiration, with control manometer, positioned on the front of the aspirator, superior to 850 mbar. Easy to carry thanks to its particular and compact dimension, the anatomic handle and its limited weight of 4,6 Kg. New case design, made of ABS/UL high resistance anti-shock with dimension of cm 34x11x24. It is supplied with 1000 ml disposable liner, green MeBer bag, 12V battery charger, suction rigid probe with hose and rechargeable bracket 10G.

The "CORMOS" aspirator for secretions is manufactured according to the ISO 10079-1:2009 product normative, IIb class, EN 1789 and is certified CE 0123, according to the directive 93/42/CEE for medical devices and european directive 104 (automotive).

### MAIN FEATURES

Width (cm)	34
Height (cm)	24
Depth (cm.)	11
Weight (kg)	4,6

### SEE ALSO



**Cormos 11010/FA**  
Portable aspirator for secretions  
with autoclavable canister

---

## ACCESSORIES

---



**11025**

Serres disposable liners 1000 ml (x  
12 pz./pcs.)



**11024**

Serres disposable liners 1000 ml (x 6  
pz./pcs.)



**11034**

Battery charger for Cormos 2 poles



**11023**

10g certified wall bracket with  
integrated charging system for  
CORMOS and CORMOS PLUS

---

## STANDARD ACCESSORIES

---



**11014**

Suction rigid probe and hose for  
"CORMOS" and "CORMOS PLUS"  
aspirators



**11016**

Green carrying bag for suction unit